

Org: pratyakshhhhhhhhhh

Phone No: 5684768946

ECG Report

DECK MOUNT

Patient Details

Name:	HARSH THAKURR	Age:	25	Gender:	Other
Date:	2025-12-18	Time:	10:24:27		

Report Overview

Maximum Heart Rate:	100 bpm
Minimum Heart Rate:	100 bpm
Average Heart Rate:	100 bpm

OBSERVATION

Interval Names	Observed Values	Standard Range
Heart Rate	100 bpm	60-100
PR Interval	146 ms	120 ms - 200 ms
QRS Complex	58 ms	70 ms - 120 ms
QRS Axis	--	Normal
QT Interval	272 ms	300 ms - 450 ms
QTcB (Bazett)	292 ms	300 ms - 450 ms
QTcF (Fridericia)	322 ms	300 ms - 450 ms
ST Deviation (J+60 ms)	--	Report in mV

ECG Report Conclusion

S.No.	Conclusion
1	Measured Values: HR 100 bpm; PR 146 ms; QRS 58 ms; QT 272 ms; QTc 292 ms; Axis --
2	Normal sinus rhythm (HR \approx 100 bpm)
3	QTcB (Bazett): 292 ms (Normal); QTcF (Fridericia): 322 ms (0.322 s)
4	ST deviation: 0.00 mV (J+60ms); report only as deviation
5	Automated interpretation (conservative): Normal unless measurements suggest otherwise
6	Acquisition: Sampling rate 500 Hz; Gain 10.0 mm/mV; Paper speed 25.0 mm/s
7	This is an automated ECG analysis and must be reviewed by a qualified physician.

Name: HARSH THAKURR

HR : 100 bpm

P/QRS/T : 132/29/54°

Age: 25

PR : 146 ms

RV5/SV1 : 1.326 mV/-0.658 mV

Gender: Other

QRS : 58 ms

RV5+SV1 : 1.984 mV

RR : 600 ms

QTCF : 322 ms

QT : 272 ms

ST : 0 ms

QTc : 292 ms

25.0 mm/s 0.5~

ANSWER

100 100 100 100 100 100 100 100 100 100

DECK MOUNT

Date: 2025-12-18

Time: 10:24:27

The figure displays twelve ECG leads (I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6) arranged vertically. Each lead shows a regular sinus rhythm with a rate of approximately 80 bpm. The strips are plotted on red grid paper. Lead I shows a narrow QRS complex with a small r wave. Lead II shows a narrow QRS complex with a prominent r wave. Lead III shows a narrow QRS complex with a large r wave. Lead aVR shows a narrow QRS complex with a small r wave. Lead aVL shows a narrow QRS complex with a small r wave. Lead aVF shows a narrow QRS complex with a small r wave. Lead V1 shows a narrow QRS complex with a small r wave. Lead V2 shows a narrow QRS complex with a small r wave. Lead V3 shows a narrow QRS complex with a small r wave. Lead V4 shows a narrow QRS complex with a small r wave. Lead V5 shows a narrow QRS complex with a small r wave. Lead V6 shows a narrow QRS complex with a small r wave.

◆ CONCLUSION ◆

1. Measured Values: HR 100 bpm; P... 2. Normal sinus rhythm (HR \approx 100 ...
3. QTcB (Bazett): 292 ms (Normal)... 4. ST deviation: 0.00 mV (J+60ms)...
5. Automated interpretation (cons... 6. Acquisition: Sampling rate 500...
7. This is an automated ECG analy...